



# **AMERICAN INTERNATIONAL UNIVERSITY-BANGLADESH**

## **CHEMISTRY-LAB SECTION: M**

**DIRECTED BY:**

**PRESENTED BY: GROUP-7**

<b>NAME</b>	<b>ID</b>
<b>SAYEDA TANJILA AHMED (APURBA)</b>	<b>*****</b>

## EXPERIMENT- 6


### **Estimation of copper (Cu) contained in a supplied solution of copper salt by iodometric method**

#### INTRODUCTION


The amount of copper ions in a supplied solution of copper salt will be measured by iodometric oxidation-reduction titration.  $\text{K}_2\text{Cr}_2\text{O}_7$  is a primary standard substance, an oxidizing agent and  $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$  is a secondary standard substance, a reducing agent.

#### REACTION

The balanced reactions are as follows:



#### **Redox Half Reactions**



## **APPARATUS**

- **BURETTE (50ML),**
- **PIPETTE (10ML),**
- **CONICAL FLASK (250ML),**
- **VOLUMETRIC FLASK (100ML),**
- **WATCH GLASS,**
- **PIPETTE FILLER,**
- **DROPPER,**
- **STAND AND CLAMP.**





## **REQUIRED CHEMICALS**


- **12% KI (Potassium iodide) solution,**
- **Copper salt solution**
- **$\text{NaHCO}_3$  (Sodium bicarbonate) solid,**
- **6M  $\text{NH}_4\text{OH}$  (Ammonium hydroxide )**
- **Conc. HCl (Hydrochloric acid),**
- **Conc.  $\text{CH}_3\text{COOH}$  (Acetic acid)**
- **Standard  $\text{K}_2\text{Cr}_2\text{O}_7$  (Potassium dichromate) solution,**
- **10%  $\text{NH}_4\text{SCN}$  (Ammonium thiocyanate) solution,**
- **$\text{Na}_2\text{S}_2\text{O}_3$  (Sodium thiosulphate) solution,**
- **Starch solution**







## PROCEDURE



**Figure: Estimation of copper (Cu) contained in a supplied solution of copper salt by iodometric method**

**Reference: (Modified version) created using ChatGPT (OpenAI,2025)**

## **EXPERIMENTAL DATA:**

**Table-1: Standardization of supplied  $\text{Na}_2\text{S}_2\text{O}_3$  solution against standard  $\text{K}_2\text{Cr}_2\text{O}_7$  solution by oxidation-reduction titration.**

No. of reading	Vol. of $\text{K}_2\text{Cr}_2\text{O}_7$ (in mL)	Vol. of $\text{Na}_2\text{S}_2\text{O}_3$ (burette reading) (in mL)			Mean (in mL)
		Initial	Final	Difference	
1	10	0	10.9	10.9	10.3
2	10	10.9	20.6	9.7	


**Table-2: Determination of the amount of copper in a supplied solution of blue vitriol by iodometric method.**

No. of reading	Vol. of Copper salt solution (in mL)	Vol. of $\text{Na}_2\text{S}_2\text{O}_3$ (burette reading) (in mL)			Mean (in mL) (V)
		Initial	Final	Difference	
1	10	0	9.1	9.1	8.73
2	10	9.1	17.7	8.6	
3	10	17.7	26.2	8.5	

## **CALCULATIONS:**

$$\text{Strength of } K_2Cr_2O_7 \text{ solution} = \frac{\text{Weight taken (in gm)} \times 0.1}{0.49} \text{ (N)} = \frac{0.5 \times 0.1}{0.49} = 0.102 \text{ N}$$

$$\begin{aligned} \text{Strength of supplied } Na_2S_2O_3 \text{ solution (S): } V_{\text{thio}} \times N_{\text{thio}} &= V_{\text{dichromate}} \times N_{\text{dichromate}} \\ \text{or, } N_{\text{thio}} &= \frac{10 \times 0.102}{10.3} = 0.099 \text{ N} \end{aligned}$$



$$\begin{aligned} \text{Amount of copper ions in 10 mL of copper salt solution} &= 0.06354 \times V \times S \text{ gm} \\ &= 0.06354 \times 8.73 \times 0.099 \text{ gm} \\ &= 0.054 \text{ gm.} \end{aligned}$$

$$\begin{aligned} \text{Amount of copper ions in 500 mL of copper salt solution} &= 0.06354 \times V \times S \times 50 \text{ gm} \\ &= 0.06354 \times 8.73 \times 0.099 \times 50 \text{ gm} \\ &= 2.70 \text{ gm.} \end{aligned}$$

$$\text{Observe value of Cu}^{2+} \text{ (in 500mL solution)} = 2.70 \text{ gm.}$$

$$\text{Known value of Cu}^{2+} \text{ (in 500mL solution)} = 2.543 \text{ gm.}$$





## **PERCENTAGE OF ERROR:**

$$\left| \frac{\text{Known value} - \text{Observed value}}{\text{Known value}} \right| \times 100\%$$

$$= \left| \frac{2.543 - 2.70}{2.543} \right| \times 100\%$$

$$= 6.17\%$$

## **RESULTS:**

***Amount of copper ions in 500 mL of copper salt solution = 2.70 gm***



# DISCUSSIONS

## CHEMISTRY LAB PRECAUTIONS



Wear gloves



Know the hazards



Wear protective eyewear



Handle chemicals properly

## ***POSSIBLE ERRORS***



In this experiment, the percentage of error obtained was 6.17%, which may be occurred due to lack of precautions or possible errors. The precautions taken for this experiment and possible errors that could have been occurred are as follows:-

The background of the slide features a blurred laboratory setting. On the left, there is a large Erlenmeyer flask containing a yellow liquid. In the center, a small vial with a cork stopper is visible. To the right, a beaker is partially filled with a red liquid. The overall scene is softly lit, with some light flares and bubbles visible in the upper corners.

### **PRECAUTIONS TAKEN:**

- **The burette was checked for air bubbles in the nozzle before starting the titration to ensure correct volume delivery.**
- **The flask was gently swirled continuously during titration to ensure uniform mixing of reactants.**
- **All glassware such as burette, pipette, and conical flask were properly cleaned with distilled water before use to avoid wrong volume.**
- **The sodium thiosulphate solution was protected from air and light as it can decompose on exposure, which may affect its strength.**
- **The titration was carried out slowly near the end point to ensure accurate detection of the color change.**





### **POSSIBLE ERRORS:**

- **Impurities present in the supplied copper salt solution may interfere with the liberation of iodine.**
- **Inaccurate standardization of sodium thiosulphate solution directly affects the calculated copper content.**
- **Parallax error may occur while reading the burette, which can cause slight variation in volume measurement.**
- **Incomplete mixing of the solution during titration may result in uneven reaction and inaccurate end point detection.**
- **Excess addition of ammonium hydroxide or acetic acid may disturb the reaction conditions, affecting the final result.**





**AMERICAN INTERNATIONAL  
UNIVERSITY-BANGLADESH**

**ON BEHALF OF  
GROUP-7**

**THANK YOU**

